

STIP1 antibody

Product Information

Catalog No.:	FNab08323
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Mediates the association of the molecular chaperones HSC70 and HSP90(HSPCA and HSPCB).

Immunogen information

Immunogen:	stress-induced-phosphoprotein 1
Synonyms:	Stress-induced-phosphoprotein 1 (STI1) Hsc70/Hsp90-organizing protein (Hop) Renal carcinoma antigen NY-REN-11 Transformation-sensitive protein IEF SSP 3521 STIP1
Observed MW:	61 kDa
Uniprot ID :	P31948

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab08323(STIP1 antibody) at dilution of 1:800