

STK31 antibody

Product Information

Catalog No.:	FNab08330
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	serine/threonine kinase 31
Synonyms:	Serine/threonine-protein kinase 31 Serine/threonine-protein kinase NYD-SPK Sugen kinase 396 (SgK396) STK31 SGK396
Observed MW:	115 kDa
Uniprot ID :	Q9BXU1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-STK31 (IP:FNab08330, 4ug;
Detection:FNab08330 1:1000) with mouse testis
tissue lysate 2400ug.



mouse testis tissue were subjected to SDS PAGE
followed by western blot with FNab08330(STK31
Antibody) at dilution of 1:600