

STK32A antibody

Product Information

Catalog No.:	FNab08331
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	serine/threonine kinase 32A
Synonyms:	Serine/threonine-protein kinase 32A Yet another novel kinase 1 STK32A YANK1
Observed MW:	46 kDa
Uniprot ID :	Q8WU08

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab08331(STK32A antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08331(STK32A antibody) at dilution of 1:800