

## STK32B antibody

### Product Information

Catalog No.:	FNab08332
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	serine/threonine kinase 32B
Synonyms:	Serine/threonine-protein kinase 32B Yet another novel kinase 2 STK32B YANK2
Observed MW:	45-48 kDa
Uniprot ID :	Q9NY57

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab08332(STK32B antibody) at dilution of 1:500