

STK32C antibody

Product Information

Catalog No.:	FNab08333
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	serine/threonine kinase 32C
Synonyms:	Serine/threonine-protein kinase 32C PKE Yet another novel kinase 3 STK32C YANK3
Observed MW:	60 kDa
Uniprot ID :	Q86UX6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human heart tissue were subjected to SDS PAGE followed by western blot with FNab08333(STK32C antibody) at dilution of 1:500