

STK38L antibody

Product Information

Catalog No.:	FNab08338
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the regulation of structural processes in differentiating and mature neuronal cells.

Immunogen information


Immunogen:	serine/threonine kinase 38 like
Synonyms:	Serine/threonine-protein kinase 38-like NDR2 protein kinase Nuclear Dbf2-related kinase 2 STK38L KIAA0965 NDR2
Observed MW:	54 kDa
Uniprot ID :	Q9Y2H1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded rat heart using FNab08338(STK38L antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab08338(STK38L antibody) at dilution of 1:800