

STOML1 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab08345 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information


| | |
|--------------|--|
| Immunogen: | stomatin(EPB72)-like 1 |
| Synonyms: | Stomatin-like protein 1 (SLP-1) EPB72-like protein 1 Protein unc-24 homolog Stomatin-related protein (STORP) STOML1 SLP1 UNC24 |
| Observed MW: | 45 kDa |
| Uniprot ID : | Q9UBI4 |

Application

| | |
|-----------------------|-----------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:20-1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded human gliomas using FNab08345(STOML1 antibody) at dilution of 1:100



A375 cells were subjected to SDS PAGE followed by western blot with FNab08345(STOML1 antibody) at dilution of 1:300