

STOML3 antibody

Product Information

Catalog No.:	FNab08348
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for the function of many mechanoreceptors. Modulate mechanotransduction channels and acid-sensing ion channels(ASIC) proteins. Potentiates PIEZO1 and PIEZO2 function by increasing their sensitivity to mechanical stimulations.

Immunogen information

Immunogen:	stomatin(EPB72)-like 3
Synonyms:	Stomatin-like protein 3 (SLP-3) STOML3
Observed MW:	30 kDa
Uniprot ID :	Q8TAV4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08348(STOML3 antibody) at dilution of 1:500