

STRBP antibody

Product Information

Catalog No.:	FNab08352
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in spermatogenesis and sperm function. Plays a role in regulation of cell growth. Binds to double-stranded DNA and RNA. Binds most efficiently to poly(I:C) RNA than to poly(dI:dC) DNA. Binds also to single-stranded poly(G) RNA. Binds non-specifically to the mRNA PRM1 3'-UTR and adenovirus VA RNA(By similarity).

Immunogen information

Immunogen:	spermatid perinuclear RNA binding protein
Synonyms:	Spermatid perinuclear RNA-binding protein STRBP SPNR
Observed MW:	74 kDa
Uniprot ID :	Q96SI9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab08352(STRBP antibody) at dilution of 1:500