

STS antibody

Product Information

Catalog No.:	FNab08354
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Conversion of sulfated steroid precursors to estrogens during pregnancy.

Immunogen information


Immunogen:	steroid sulfatase(microsomal), isozyme S
Synonyms:	Steryl-sulfatase Arylsulfatase C (ASC) Estrone sulfatase Steroid sulfatase Steryl-sulfate sulfohydrolase STS ARSC1
Observed MW:	63 kDa-65 kDa
Uniprot ID :	P08842

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNab08354(STS antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab08354(STS antibody) at dilution of 1:500