

## STT3B antibody

### Product Information

Catalog No.:	FNab08356
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	STT3, subunit of the oligosaccharyltransferase complex, homolog B
Synonyms:	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit STT3B (Oligosaccharyl transferase subunit STT3B, STT3-B) Source of immunodominant MHC-associated peptides homolog STT3B SIMP
Observed MW:	94 kDa
Uniprot ID :	Q8TCJ2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab08356( STT3B Antibody) at dilution of 1:200



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08356(STT3B antibody) at dilution of 1:2000