

## SYT3 antibody

### Product Information

Catalog No.:	FNab08429
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	synaptotagmin III
Synonyms:	Synaptotagmin-3 Synaptotagmin III (SytIII) SYT3
Observed MW:	63 kDa
Uniprot ID :	Q9BQG1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab08429(Synaptotagmin-3 antibody) at dilution of 1:300