

## SYNJ1 antibody

### Product Information

Catalog No.:	FNab08438
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Inositol 5-phosphatase which has a role in clathrin-mediated endocytosis.

### Immunogen information

Immunogen:	synaptojanin 1
Synonyms:	Synaptojanin-1 Synaptic inositol 1,4,5-trisphosphate 5-phosphatase 1 SYNJ1 KIAA0910
Observed MW:	140 kDa
Uniprot ID :	O43426

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab08438( SYNJ1 Antibody) at dilution of 1:4000