

STX18 antibody

Product Information

Catalog No.:	FNab08449
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	syntaxin 18
Synonyms:	Syntaxin-18 Cell growth-inhibiting gene 9 protein STX18
Observed MW:	42 kDa
Uniprot ID :	Q9P2W9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2400; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human endometrial cancer using FNab08449(STX18 antibody) at dilution of 1:100



human liver tissue were subjected to SDS PAGE followed by western blot with FNab08449(STX18 antibody) at dilution of 1:500