

## STX4 antibody

### Product Information

Catalog No.:	FNab08452
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	syntaxin 4
Synonyms:	Syntaxin-4 Renal carcinoma antigen NY-REN-31 STX4 STX4A
Observed MW:	37 kDa
Uniprot ID :	Q12846

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A375 cells were subjected to SDS PAGE followed by western blot with FNab08452( Syntaxin 4 antibody) at dilution of 1:1000