

## STX8 antibody

### Product Information

Catalog No.:	FNab08456
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	syntaxin 8
Synonyms:	Syntaxin-8 STX8
Observed MW:	27 kDa
Uniprot ID :	Q9UNK0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, WB, FC
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000;IF: 1:10-1:100
Image:	



Immunofluorescent analysis of HepG2 cells, using STX8 antibody FNab08456 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



human heart tissue were subjected to SDS PAGE followed by western blot with FNab08456(Syntaxin 8 antibody) at dilution of 1:500