

SYTL3 antibody

Product Information

Catalog No.:	FNab08461
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May act as Rab effector protein and play a role in vesicle trafficking. Binds phospholipids in the presence of calcium ions(By similarity).

Immunogen information

Immunogen:	synaptotagmin-like 3
Synonyms:	Synaptotagmin-like protein 3 Exophilin-6 SYTL3 SLP3
Observed MW:	69 kDa
Uniprot ID :	Q4VX76

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab08461(SYTL3 antibody) at dilution of 1:500