

TACO1 antibody

Product Information

Catalog No.:	FNab08470
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	coiled-coil domain containing 44
Synonyms:	Translational activator of cytochrome c oxidase 1 Coiled-coil domain-containing protein 44 Translational activator of mitochondrially-encoded cytochrome c oxidase I TACO1 CCDC44
Observed MW:	30-32 kDa
Uniprot ID :	Q9BSH4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab08470(CCDC44 antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab08470(CCDC44 Antibody) at dilution of 1:600