

## TACR2 antibody

### Product Information

Catalog No.:	FNab08471
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


This is a receptor for the tachykinin neuropeptide substance K(neurokinin A). It is associated with G proteins that activate a phosphatidylinositol-calcium second messenger system. The rank order of affinity of this receptor to tachykinins is: substance K > neuromedin-K > substance P.

### Immunogen information

Immunogen:	tachykinin receptor 2
Synonyms:	Substance-K receptor (SKR) NK-2 receptor (NK-2R) Neurokinin A receptor Tachykinin receptor 2 TACR2 NK2R NKNAR TAC2R
Observed MW:	44 kDa
Uniprot ID :	P21452

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08471( TACR2 Antibody) at dilution of 1:600