

## TADA1 antibody

### Product Information

Catalog No.:	FNab08473
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	transcriptional adaptor 1(HFI1 homolog, yeast)-like
Synonyms:	Transcriptional adapter 1 SPT3-associated factor 42 (STAF42) Transcriptional adapter 1-like protein TADA1 TADA1L
Observed MW:	38 kDa
Uniprot ID :	Q96BN2

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab08473(TADA1L Antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08473(TADA1L Antibody) at dilution of 1:300