

TAF13 antibody

Product Information

Catalog No.:	FNab08479
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


TFIID beta-specific TAFII.

Immunogen information

Immunogen:	TAF13 RNA polymerase II, TATA box binding protein(TBP)-associated factor, 18kDa
Synonyms:	Transcription initiation factor TFIID subunit 13 Transcription initiation factor TFIID 18 kDa subunit (TAF(II)18, TAFII-18, TAFII18) TAF13 TAF2K TAFII18
Observed MW:	14-18 kDa
Uniprot ID :	Q15543

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08479(TAF13 Antibody) at dilution of 1:600