

## TAF1-Specific antibody

### Product Information

Catalog No.:	FNab08484
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	TAF1 RNA polymerase II, TATA box binding protein(TBP)-associated factor, 250kDa
Synonyms:	Transcription initiation factor TFIID subunit 1 Cell cycle gene 1 protein TBP-associated factor 250 kDa (p250) Transcription initiation factor TFIID 250 kDa subunit (TAF(II)250, TAFII-250, TAFII250) TAF1 BA2R CCG1 CCGS TAF2A
Observed MW:	50 kDa
Uniprot ID :	P21675

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab08484(TAF1-Specific antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab08484(TAF1-Specific antibody) at dilution of 1:300