

TAF6L antibody

Product Information

Catalog No.:	FNab08487
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Functions as a component of the PCAF complex. The PCAF complex is capable of efficiently acetylating histones in a nucleosomal context. The PCAF complex could be considered as the human version of the yeast SAGA complex.

Immunogen information

Immunogen:	TAF6-like RNA polymerase II, p300/CBP-associated factor(PCAF)-associated factor, 65kDa
Synonyms:	TAF6-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 6L (TAF6L) PCAF-associated factor 65-alpha (PAF65-alpha) TAF6L PAF65A
Observed MW:	68 kDa
Uniprot ID :	Q9Y6J9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab08487(TAF6L antibody) at dilution of 1:500