

TLN1 antibody

Product Information

Catalog No.:	FNab08495
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	talin 1
Synonyms:	Talin-1 TLN1 KIAA1027 TLN
Observed MW:	270 kDa, 230 kDa
Uniprot ID :	Q9Y490

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab08495(TLN1 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08495(TLN1 antibody) at dilution of 1:500