

## TRIOBP antibody

### Product Information

Catalog No.:	FNab08500
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	TRIO and F-actin binding protein
Synonyms:	TRIO and F-actin-binding protein Protein Tara TRF1-associated protein of 68 kDa Trio-associated repeat on actin TRIOBP KIAA1662 TARA
Observed MW:	68-75 kDa
Uniprot ID :	Q9H2D6

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab08500(TARA antibody) at dilution of 1:1000