

TBC1D1 antibody

Product Information

Catalog No.:	FNab08509
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May act as a GTPase-activating protein for Rab family protein(s). May play a role in the cell cycle and differentiation of various tissues. Involved in the trafficking and translocation of GLUT4-containing vesicles and insulin-stimulated glucose uptake into cells(By similarity).

Immunogen information


Immunogen:	TBC1(tre-2/USP6, BUB2, cdc16) domain family, member 1
Synonyms:	TBC1 domain family member 1 TBC1D1 KIAA1108
Observed MW:	150 kDa
Uniprot ID :	Q86TI0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney slide using FNab08509(TBC1D1 Antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab08509(TBC1D1 antibody) at dilution of 1:1000