

## TBCB antibody

### Product Information

Catalog No.:	FNab08517
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds to alpha-tubulin folding intermediates after their interaction with cytosolic chaperonin in the pathway leading from newly synthesized tubulin to properly folded heterodimer. Involved in regulation of tubulin heterodimer dissociation. May function as a negative regulator of axonal growth.

### Immunogen information


Immunogen:	tubulin folding cofactor B
Synonyms:	Tubulin-folding cofactor B Cytoskeleton-associated protein 1 Cytoskeleton-associated protein CKAP1 Tubulin-specific chaperone B TBCB CG22 CKAP1
Observed MW:	32 kDa
Uniprot ID :	Q99426

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-TBCB (IP:FNab08517, 3ug;  
Detection:FNab08517 1:500) with HeLa cells  
lysate 1200ug.



HeLa cells were subjected to SDS PAGE followed  
by western blot with FNab08517(TBCB antibody)  
at dilution of 1:500