

TBCEL antibody

Product Information

Catalog No.:	FNab08521
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a regulator of tubulin stability.

Immunogen information

Immunogen:	tubulin folding cofactor E-like
Synonyms:	Tubulin-specific chaperone cofactor E-like protein (EL) Leucine-rich repeat-containing protein 35 TBCEL LRRC35
Observed MW:	40-48 kDa
Uniprot ID :	Q5QJ74

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08521(TBCEL Antibody) at dilution of 1:3000