

TBRG1 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab08531 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


Acts as a growth inhibitor. Can activate p53/TP53, causes G1 arrest and collaborates with CDKN2A to restrict proliferation, but does not require either protein to inhibit DNA synthesis. Redistributes CDKN2A into the nucleoplasm. Involved in maintaining chromosomal stability.

Immunogen information

| | |
|--------------|---|
| Immunogen: | transforming growth factor beta regulator 1 |
| Synonyms: | Transforming growth factor beta regulator 1 Nuclear interactor of ARF and Mdm2 TBRG1 NIAM |
| Observed MW: | 45 kDa |
| Uniprot ID : | Q3YBR2 |

Application

| | |
|--------------------------|--------------|
| Reactivity: | Human |
| Tested Application: | ELISA, WB |
| Recommended dilution:WB: | 1:200-1:2000 |
| Image: | |



293T cells were subjected to SDS PAGE followed by western blot with FNab08531(TBRG1 antibody) at dilution of 1:1000