

## TBX15 antibody

### Product Information

Catalog No.:	FNab08534
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probable transcriptional regulator involved in the development of the skeleton of the limb, vertebral column and head. Acts by controlling the number of mesenchymal precursor cells and chondrocytes(By similarity).

### Immunogen information

Immunogen:	T-box 15
Synonyms:	T-box transcription factor TBX15 (T-box protein 15) T-box transcription factor TBX14 (T-box protein 14) TBX15 TBX14
Observed MW:	55 kDa
Uniprot ID :	Q96SF7

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08534(TBX15 Antibody) at dilution of 1:600