

## TCF12 antibody

### Product Information

Catalog No.:	FNab08552
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	transcription factor 12
Synonyms:	Transcription factor 12 (TCF-12) Class B basic helix-loop-helix protein 20 (bHLHb20) DNA-binding protein HTF4 E-box-binding protein Transcription factor HTF-4 TCF12 BHLHB20 HEB HTF4
Observed MW:	73 kDa
Uniprot ID :	Q99081

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab08552(TCF12 antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab08552(TCF12 antibody) at dilution of 1:300