

TCN2 antibody

Product Information

Catalog No.:	FNab08561
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Primary vitamin B12-binding and transport protein. Delivers cobalamin to cells.

Immunogen information


Immunogen:	transcobalamin II; macrocytic anemia
Synonyms:	Transcobalamin-2 (TC-2) Transcobalamin II (TC II, TCII) TCN2 TC2
Observed MW:	48 kDa
Uniprot ID :	P20062

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse small intestine using FNa08561(TCN2 antibody) at dilution of 1:50



A375 cells were subjected to SDS PAGE followed by western blot with FNa08561(TCN2 antibody) at dilution of 1:100