

TCP10L antibody

Product Information

Catalog No.:	FNab08565
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation. Has in vitro transcription inhibition activity. Acts as a tumor suppressor in hepatocellular carcinoma(HCC) cells.

Immunogen information


Immunogen:	t-complex 10(Mouse)-like
Synonyms:	T-complex protein 10A homolog 1 (T-complex protein 10A-1, TCP10A-1) TCP10-like TCP10L
Observed MW:	24 kDa
Uniprot ID :	Q8TDR4

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab08565(TCP10L antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08565(TCP10L antibody) at dilution of 1:200