

TCTN3 antibody

Product Information

Catalog No.:	FNab08569
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Part of the tectonic-like complex which is required for tissue-specific ciliogenesis and may regulate ciliary membrane composition(By similarity). May be involved in apoptosis regulation. Necessary for signal transduction through the sonic hedgehog(Shh) signaling pathway.

Immunogen information


Immunogen:	tectonic family member 3
Synonyms:	Tectonic-3 TCTN3 C10orf61 TECT3
Observed MW:	66 kDa
Uniprot ID :	Q6NUS6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung using FNab08569(TCTN3 antibody) at dilution of 1:50



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08569(TCTN3 Antibody) at dilution of 1:1000