

## TNR antibody

### Product Information

Catalog No.:	FNab08595
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	tenascin R(restrictin, janusin)
Synonyms:	Tenascin-R (TN-R) Janusin Restrictin TNR
Observed MW:	160 kDa
Uniprot ID :	Q92752

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab08595(TNR-Specific antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08595(TNR-Specific antibody) at dilution of 1:300