

TEX101 antibody

Product Information

Catalog No.:	FNab08605
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in signal transduction and promote protein tyrosine phosphorylation.

Immunogen information

Immunogen:	testis expressed 101
Synonyms:	Testis-expressed protein 101 Cell surface receptor NYD-SP8 Scleroderma-associated autoantigen Spermatogenesis-related gene protein TEX101 SGRG
Observed MW:	None
Uniprot ID :	Q9BY14

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab08605 (TEX101 antibody) at dilution of 1:50