

## TEX13A antibody

### Product Information

Catalog No.:	FNab08608
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	testis expressed 13A
Synonyms:	Testis-expressed protein 13A TEX13A
Observed MW:	46 kDa
Uniprot ID :	Q9BXU3

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab08608(TEX13A antibody) at dilution of 1:100



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab08608(TEX13A antibody) at dilution of 1:500