

TFAP2E antibody

Product Information

Catalog No.:	FNab08620
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Sequence-specific DNA-binding protein that interacts with inducible viral and cellular enhancer elements to regulate transcription of selected genes. AP-2 factors bind to the consensus sequence 5'-GCCNNNGGC-3' and activate genes involved in a large spectrum of important biological functions including proper eye, face, body wall, limb and neural tube development. They also suppress a number of genes including MCAM/MUC18, C/EBP alpha and MYC. AP-2-epsilon may play a role in the development of the CNS and in cartilage differentiation(By similarity).


Immunogen information

Immunogen:	transcription factor AP-2 epsilon(activating enhancer binding protein 2 epsilon)
Synonyms:	Transcription factor AP-2-epsilon (AP2-epsilon) Activating enhancer-binding protein 2-epsilon TFAP2E
Observed MW:	46 kDa
Uniprot ID :	Q6VUC0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000

Image:



A375 cells were subjected to SDS PAGE followed by western blot with FNab08620(TFAP2E Antibody) at dilution of 1:600