

GTF2I antibody

Product Information

Catalog No.:	FNab08630
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	general transcription factor II, i
Synonyms:	General transcription factor II-I (GTFII-I, TFII-I) Bruton tyrosine kinase-associated protein 135 (BAP-135, BTK-associated protein 135) SRF-Phox1-interacting protein (SPIN) Williams-Beuren syndrome chromosomal region 6 protein GTF2I BAP135 WBSCR6
Observed MW:	135 kDa
Uniprot ID :	P78347

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, FC, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab08630(GTF2I Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08630(GTF2I antibody) at dilution of 1:800