

TGDS antibody

Product Information

Catalog No.:	FNab08636
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	TDP-glucose 4,6-dehydratase
Synonyms:	dTDP-D-glucose 4,6-dehydratase TGDS
Observed MW:	45-48 kDa
Uniprot ID :	O95455

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human bladder using FNab08636(TGDS antibody) at dilution of 1:50



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab08636(TGDS antibody) at dilution of 1:300