

TGM4 antibody

Product Information

Catalog No.:	FNab08651
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Associated with the mammalian reproductive process. Catalyzes the cross-linking of proteins and the conjugation of polyamines to specific proteins in the seminal tract.

Immunogen information


Immunogen:	transglutaminase 4(prostate)
Synonyms:	Protein-glutamine gamma-glutamyltransferase 4 Fibrinoligase Prostate transglutaminase Prostate-specific transglutaminase Transglutaminase P (TG(P), TGP, TGase P) Transglutaminase-4 (TGase-4) TGM4
Observed MW:	68-70 kDa
Uniprot ID :	P49221

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab08651(TGM4 antibody) at dilution of 1:100



human testis tissue were subjected to SDS PAGE followed by western blot with FNab08651(TGM4 antibody) at dilution of 1:300