

TGS1 antibody

Product Information

Catalog No.:	FNab08654
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information


Immunogen:	trimethylguanosine synthase homolog
Synonyms:	Trimethylguanosine synthase CLL-associated antigen KW-2 Cap-specific guanine-N2 methyltransferase Hepatocellular carcinoma-associated antigen 137 Nuclear receptor coactivator 6-interacting protein PRIP-interacting protein with methyltransferase motif (PIMT, PIPMT) TGS1 HCA137 NCOA6IP PIMT
Observed MW:	97 kDa
Uniprot ID :	Q96RS0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemical of paraffin-embedded human kidney using FNaB08654(TSG1 antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNaB08654(TSG1 antibody) at dilution of 1:1000