

THAP2 antibody

Product Information

Catalog No.:	FNab08661
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	THAP domain containing, apoptosis associated protein 2
Synonyms:	THAP domain-containing protein 2 THAP2
Observed MW:	26-30 kDa
Uniprot ID :	Q9H0W7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08661(THAP2 antibody) at dilution of 1:600