

ACOT13 antibody

Product Information

Catalog No.:	FNab08662
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	thioesterase superfamily member 2
Synonyms:	Acyl-coenzyme A thioesterase 13 (Acyl-CoA thioesterase 13) Hotdog-fold thioesterase superfamily member 2 Palmitoyl-CoA hydrolase Thioesterase superfamily member 2 (THEM2) Acyl-coenzyme A thioesterase 13, N-terminally processed ACOT13 THEM2
Observed MW:	16 kDa
Uniprot ID :	Q9NPJ3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab08662(THEM2 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08662(THEM2 antibody) at dilution of 1:500