

## THNSL2 antibody

### Product Information

Catalog No.:	FNab08665
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Isoform 1: Acts as a catabolic phospho-lyase on both gamma-and beta-phosphorylated substrates. Degrades O-phospho-threonine(PThr) to alpha-ketobutyrate, ammonia and phosphate(By similarity). Isoform SOFAT: Potent inducer of osteoblastic production of IL6. May act to exacerbate inflammation and/or bone turnover under inflammatory conditions.

### Immunogen information

Immunogen:	threonine synthase-like 2
Synonyms:	Threonine synthase-like 2 (TSH2) Secreted osteoclastogenic factor of activated T-cells (SOFAT) THNSL2
Observed MW:	55-60 kDa
Uniprot ID :	Q86YJ6

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08665(THNSL2 antibody) at dilution of 1:1000