

THOC1 antibody

Product Information

Catalog No.:	FNab08666
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


HPR1 is part of the TREX (transcription/export) complex, which includes TEX1 (MIM 606929), THO2 (MIM 300395), ALY (MIM 604171), and UAP56 (MIM 142560).

Immunogen information


Immunogen:	THO complex 1
Synonyms:	THO complex subunit 1 (Tho1) Nuclear matrix protein p84 (p84N5) hTREX84 THOC1 HPR1
Observed MW:	80 kDa
Uniprot ID :	Q96FV9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab08666(THOC1 antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08666(THOC1 antibody) at dilution of 1:1000