

## THUMPD1 antibody

### Product Information

Catalog No.:	FNab08679
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Functions as a tRNA-binding adapter to mediate NAT10-dependent tRNA acetylation(PubMed:25653167).

### Immunogen information


Immunogen:	THUMP domain containing 1
Synonyms:	THUMP domain-containing protein 1 THUMPD1
Observed MW:	49 kDa
Uniprot ID :	Q9NXG2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab08679(THUMPD1 antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab08679(THUMPD1 antibody) at dilution of 1:1000