

TIMM17A antibody

Product Information

Catalog No.:	FNab08697
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information


Immunogen:	translocase of inner mitochondrial membrane 17 homolog A (yeast)
Synonyms:	Mitochondrial import inner membrane translocase subunit Tim17-A Inner membrane preprotein translocase Tim17a TIMM17A MIMT17 TIM17 TIM17A TIMM17
Observed MW:	18 kDa
Uniprot ID :	Q99595

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab08697(TIMM17A antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08697(TIMM17A antibody) at dilution of 1:1000