

## TIMMDC1 antibody

### Product Information

Catalog No.:	FNab08705
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Chaperone protein involved in the assembly of the mitochondrial NADH:ubiquinone oxidoreductase complex(complex I). Participates in constructing the membrane arm of complex I.

### Immunogen information

Immunogen:	chromosome 3 open reading frame 1
Synonyms:	Complex I assembly factor TIMMDC1, mitochondrial Protein M5-14 Translocase of inner mitochondrial membrane domain-containing protein 1 (TIMM domain containing-protein 1) TIMMDC1 C3orf1
Observed MW:	33 kDa
Uniprot ID :	Q9NPL8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



A375 cells were subjected to SDS PAGE followed by western blot with FNab08705(TIMMDC1 antibody) at dilution of 1:1000