

## TINAGL1 antibody

### Product Information

Catalog No.:	FNab08711
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be implicated in the adrenocortical zonation and in mechanisms for repressing the CYP11B1 gene expression in adrenocortical cells. This is a non catalytic peptidase C1 family protein(By similarity).

### Immunogen information


Immunogen:	tubulointerstitial nephritis antigen-like 1
Synonyms:	Tubulointerstitial nephritis antigen-like Glucocorticoid-inducible protein 5 Oxidized LDL-responsive gene 2 protein (OLRG-2) Tubulointerstitial nephritis antigen-related protein (TIN Ag-related protein, TIN-Ag-RP) TINAGL1 GIS5 LCN7 OLRG2 TINAGL
Observed MW:	52 kDa
Uniprot ID :	Q9GZM7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-TINAGL1 (IP:FNab08711, 4ug;  
Detection:FNab08711 1:500) with mouse heart  
tissue lysate 4000ug.



human placenta tissue were subjected to SDS  
PAGE followed by western blot with  
FNab08711(TINAGL1 antibody) at dilution of  
1:800